# Application Manual: For Use in CS6460

Fei Ding, Jonathan Hou, Jan Jendrusak, Jintong Jiang {fding33, jhou47, jjendrusak3, jjiang322}@gatech.edu

**Abstract**—This document is the manual of our application, which includes the instructions to utilize the application and explore all the features.

### 1 Overview

Our application is deployed with Heroku and you can access the application via the link: <a href="https://cs6440-dermatology-risk.herokuapp.com/">https://cs6440-dermatology-risk.herokuapp.com/</a>. With our application, patients can upload a skin image and get a score that indicates the probability of getting skin cancer. A higher score represents that the user has a higher risk of getting cancer. The user can create an account and our application will store the previous images he uploads. Also, he can use our system to submit an inquiry of sending an email to a specific dermatologist. We will introduce each feature in the following sections.

## 2 Login

There are three functions below the navigation bar: Home, Upload new images, and inquiry. To create an account and log into our application, the user need to provide a pair of username and password.




Figure 1 — Navigation Bar and the UI of login function

After creating a new account, the user is asked to provide the date of birth and sex to create a profile. This information is not required to make the model work, but it would be useful when creating an inquiry to contact a dermatologist.




Figure 2 — Profile Page

#### 3 UPLOAD NEW IMAGE

The user can choose the "upload new image" function to upload an image to diagnose. The image should be less than 200 MB and we support the image with PNG, JPG, and JPEG formats. The user could choose the date when the image is taken and also add descriptions related to the skin problems. All the images this user uploads will be stored in the database. After clicking submit, the website will show a notification that says "Success! Your image has been uploaded. Model result will show up in a few seconds," if the image is uploaded successfully. The user can go back to the home page to check the results after waiting for a few seconds. The home page will display the most recent image and match overview, which includes the predicted risk score and the risk of developing skin cancer.




Figure 3— Upload new image to run a diagnosis.




Figure 4— Match Overview and results

# 4 INQUIRY

If the user wants to contact a dermatologist to further diagnose his skin issue, he can choose the inquiry function. When creating the inquiry, our application won't share the user's profile information, including sex and age, with the dermatologist. It will send the images of the skin the user uploaded previously and the model prediction score to the dermatologist. When creating the inquiry, the user should provide the name and email address of the dermatologist. Also, the user can also provide additional notes to the doctor to help the doctor

diagnose conveniently. After clicking the "send email" button, the user can send the email to the corresponding dermatologist.




Figure 5 — Inquiry page